


Tape or
Glue


Fold

Tape or
Glue

Fold




PCG Bot™
©PCG 2020



1




2




PCG Bot Manual

- 1 color Zoovacs stick Gray and his hair Green
- 2 color the  the color of choice
- 3 color A, B, X, Y,  = Analog stick,  the color of choice
- 4 Fold the indicated  side inward
- 5 Fold the indicated side inward
- 6 Fold the indicated side inward
- 7 color in  the color of choice
- 8 cut the labeled area cut →
9. Apply Tape or Glue stick Adhesive on the indicated area

Apps



PCG Bot



System / system 1.0

PCG Bot



Vol.



Brightness

GB  200

Your BOT

